



Data sheet

Specification of the product

Article code: **1806-PH**
Description: **Wonder nebuliser, T-piece, mouthpiece with 2,10m tubing**

Statistical information:

Product class according to 93/42/EEC: II a
Customs code: 90192000
Country of origin: Poland
Packaging: each piece clean packed in a poly-bag, each 50 pieces packed in a carton-box
Box size: 580 mm x 260 mm x 270 mm
Weight net of a one packed piece: 85g
Weight of a 50 pieces box: 4,75kg

Indications for use: administer liquid medicines to the patient in the form of an aerosol by inhalation through the mouth directly into the respiratory system

Material information: Tubing: 100% medical PVC according to European Pharmacopeia 100% Latex free, 100% Phthalate free concerning PVC (production process: injection molding, extrusion);
Nebuliser top and the element inside: polypropylene
Nebuliser medicine container: polystyrene
T-piece: polypropylene
Mouthpiece: polypropylene

Labeling/Lot number: 1 label with the description in 7 languages according to the valid production procedure inside of each poly-bag;
2 labels on each carton
Each label includes the changing Lot number. The lot number includes 6 signs: the first 2 means the production year, the next the production month and the last 2 the day the begin of the production.

Expiratory date: (see Labeling)

| | |
|-----------------------------|---|
| Storage: | keep the goods in a dry and clean room acc. to the local standards for medical devices |
| Conditions of use: | temperature: +5°C to +40°C |
| Sterile: | No |
| Instruction for use: | due to easy use the product is self explanatory. The exact usage of the product depends on the type of the therapy and the instruction of a physician. This product is a single patient use product. There must be assured, the condition and the cleanliness of the product is good enough before twice use. |
| Recycling/ Wastes: | use the national standards for the destruction of PVC medical products. |
| Related documents: | ISO certificates Declaration of conformity |